



## ZOO MED'S Mail-in Rebate Card

For Zoo Med Reptibator®  
Digital Egg Incubator

### TO RECEIVE YOUR ZOO MED REBATE:

- 1) Purchase a qualifying product(s) during the promotion dates:  
01/1/18 through 12/31/19.
- 2) Remove the original UPC from the qualifying product packaging.  
Copies of the UPC do not qualify.
- 3) Mail the following items to Zoo Med:
  - Original UPC
  - A copy of your invoice with purchase(s) circled
  - This completed offer form within 30 days of the purchase date to the following address:

SEND TO: ZOO MED REBATES: 3650 Sacramento Dr. - San Luis Obispo, CA 93401

<hr/>		
First Name	Last Name	
<hr/>		
Address (sorry, we cannot deliver to P.O. boxes)		
<hr/>		
City	State	Zip/Postal Code
<hr/>		
Daytime Phone	E-mail Address	

#### TERMS AND CONDITIONS:

Your right to receive this rebate will not be earned unless you satisfy these Conditions. Failure to satisfy each of these Conditions is a rejection of this rebate. This offer cannot be combined with any other offer, unless specified in writing by Zoo Med Laboratories, Inc. This offer is available to all customers with mailing addresses in the United States only. Purchases made in or for delivery to other countries are not eligible. All rebate checks will be paid in U.S. dollars. Rebate Sponsor is not responsible for lost, late, damaged, illegible, misdirected or postage-due submissions. Your rights to this offer cannot be assigned or transferred and this offer is void where taxed, restricted or prohibited by law. Excessive submissions constitute fraud and may result in federal prosecution under the U.S. mail fraud statutes (Title 18, USC 1341 and 1342). All submitted materials become property of Zoo Med Laboratories and will not be returned. It is required that you cash the rebate check within 90 days of the issuance date on the check. Check will be mailed to qualified customers within 8 weeks of the postmark date on your request.



# \$5.00 off

Zoo Med's Reptibator®  
Digital Egg Incubator



#### FEATURES

- Digital controller with LCD display and LED heat indicator light.
- Pulse proportional thermostat for stable temperature regulation.
- Temperature control range from 59°F to 104°F (15°C to 40°C).
- Humidity display range from 10% to 95% relative humidity.
- Rigid DS with heating element.
- Programmable temperature alarm.
- Built-in memory stores settings in case of power failure.
- Clear cover allows for easy egg inspection without opening incubator.
- Multiple outer recesses molded into incubator base. Fee-rough ventilation and sponge insert allow for even evaporation and proper humidity levels within the incubator.



Heat indicator with display and LED heat indicator light.



Multiple outer recesses molded into incubator base. Fee-rough ventilation and sponge insert allow for even evaporation and proper humidity levels within the incubator.



Clear cover allows for easy egg inspection without opening incubator.



Multiple outer recesses molded into incubator base. Fee-rough ventilation and sponge insert allow for even evaporation and proper humidity levels within the incubator.

INSTRUCTIONS  
ENCLOSED